



EDUCATIONAL SERVICE
UNIT 8
SPECIAL EDUCATION
DEPARTMENT LOANABLE &
DISPOSABLE EQUIPMENT
2012-2013

LOANABLE
EQUIPMENT

Pony Walker - 2
Baby Bath Seat - 3
Corner Floor Sitter - 4
Bumbo - 5
Pilates Roller - 6
Half Round Rolls - 6
Adjustable Bench - 7
Pudgey Piglet - 8
Stacking Ring Cone - 8
Corner Seat Position Aid - 9

IF YOU ARE
INTERESTED IN
BORROWING ANY OF
OUR LOANABLE
EQUIPMENT OR
WOULD LIKE ONE OF
THE DISPOSABLE
PRODUCTS WE OFFER
PLEASE CONTACT
TONYA STROM AT
(402)887-5041 EXT. 1242
OR stromt@esu8.org

DISPOSABLE
EQUIPMENT

Grabber - 10
Scooper Plate - 11
Nosey Cup - 11
Spoons - 12
Utensil Holder - 12
Slip on Handle - 12

PONY WALKER

The Pony is a gait trainer for disabled children. The Pony has a frame shaped like an arrow with five castors. A central bar holds the chest support and the seat with pelvic support. The bar can be tilted forwards to stimulate walking. Between the legs is a plate to avoid cross stepping. The Pony is for indoor use. See website for additional information.




A

B

C

DREAM BABY BATH SEAT

The Dream Baby Premium Deluxe Bath Seat makes it much easier to bathe a child. The front support bar of the seat can easily be opened, which simplifies the positioning of baby while getting in and out of the bath seat. The seat is lightweight and durable with strong suction cups to secure the seat to the bath, plus the seat can be easily folded for convenient storage. This seat can be used for children 6 to 24 months.



RIFTON CORNER FLOOR SITTER CHAIR

This chair can be used for children whom cannot sit up on their own, and who need their shoulders and upper trunk supported. This seat provides a simple way to seat a child at floor level. In this low chair the child is independent, with hands free for play, while maintaining optimal erect trunk alignment. When the chest strap of this mobile positioning chair is loosened, the child has limited movement, while remaining well secured. This encourages trunk strength, stability, and self-righting. Sitting at a 90-degree angle, with comfortable leg extension, children can join in on the level they like best.



BUMBO BABY SEAT AND PLAY TRAY

The Bumbo Baby Seat assists babies to sit upright when babies can adequately support their own head. The Bumbo should be used on a safe ground level surface under adult supervision. The Play Tray can be attached and secured to the seat or the seat can be used on its own.



INFLATABLE PILATES ROLLER

Versatile and inflatable to desired firmness, this roller is made of durable latex-free vinyl material. The roller is excellent for core, flexibility, and balance training.



FOAM HALF-ROUND ROLLS

Lightweight foam rolls can be used for working on standing balance or with smaller children on tummy and on all fours. Two sizes available 3"D X 12"L & 3"D X 3"L.



KAYE ADJUSTABLE BENCH

This bench will provide a stable sitting surface to achieve maximum arm and hand movements for play, learning, and daily living activities. This bench is good for use at home, in school programs and for children who do not need back support in sitting but need a firm seat that can be adjusted so that the child's feet rest on the floor. This bench is also ideal for children learning to move from one level to another, such as moving from the floor to sitting, from sitting to standing, or from a seated position to their walker. The bench can be raised, lowered, or tilted quickly and easily. The bench has a user weight limit of 180 lbs.



PUDGEY THE PIGLET

Switch adapted toy, Pudgey walks, wags his tail, wrinkles his nose, and oinks to his heart's content. Good for increasing visual tracking and visual attention skills as well as cause and effect. Pudgey can be activated by a press button switch as pictured, and/or can be connected to the stacking ring cone switch pictured below. The stacking ring cone can also be used alone. When the colorful wooden rings are stacked in the proper order, the built-in music box will begin to play. This toy is great for teaching discrimination between sizes and colors and it is great for sharpening their visual-motor skills, too.



CORNER SEAT POSITIONING AID

The Pediatric Corner Seat and Positioning aid can be used anywhere – indoors, outdoors, bathtub. It is designed to promote head and upper extremity control in order to encourage the development of fine motor control. The corner seat stabilizes the shoulders, pelvis and trunk and facilitates the arms and shoulders to come forward, which increases arm movements and allows the child to work on activities in the midline. The seat features an H-Strap for anterior trunk support, a hip belt to maintain proper position and an abductor block to help maintain leg extension and reduce extensor thrust. The broad base design provides stability and helps to prevent the seat from tipping over. The seat is made of polyethylene and is easy to wash and maintain.



ARR'S GRABBER XT

As a mouthing tool, ARK's Grabber provides an excellent opportunity for the mouth to play through jaw movement, tongue movement, and oral exploration. The chewy, resilient material helps with oral defensiveness and food aversions by de-sensitizing the mouth. It is also a safe and durable alternative to tooth grinding and/or chewing on hazardous non-food items. Each Grabber is constructed of medical grade materials that are FDA compliant and contain no lead, phthalates, PVC, BPA, or latex. The resilient material is solid with no through holes or crevices that could harbor bacteria or saliva. The Grabber can be used as an alternative to biting/chewing on fingers, shirts, pencils, and other hazardous non-food items. It can also be used to decrease pica, thumb sucking, and tooth grinding, to build oral strength and tolerance, and to calm and soothe individuals through oral sensory input.



SCOOPER PLATE

The flat bottom and high rim with a reverse curve on one side help in scooping food onto a utensil without spilling over the side. The plate diameter is 6.75" and is molded of heat-resistant plastic.



NOSEY CUP

The rim of the Nosey cup has been designed so a person can drink without tipping their head back to drain the glass. It can be tipped beyond the point where a normal cup would hit the bridge of your nose. The cup is designed for people with limited neck movement, a bed patient, or in a wheelchair with a headrest. Can be microwaved for warming only. The cup we have available is 4 oz. and transparent green.



SMALL MAROON SPOON (5 PACK)

This spoon is a great for providing feeding therapy. The spoon has a one-inch bowl width and is made out of FDA approved plastic that contain NO latex, lead, BPA, or phthalates. Although the spoons are durable and strong, they must not be used with individuals who have a strong bite reflex.



KINSMAN UTENSIL HOLDER – RAINBOW

This utensil holder is an excellent eating aid for people who have trouble grasping and holding small utensils. This spring action plastic clip fits the hand snugly without pinching. Great for people with grasp limitations or trouble manipulating small utensils.



SLIP-ON UTENSIL HANDLE

The utensil holder designed for use by individuals with upper extremity or grasping disabilities or spinal cord injury. The outside surface of this holder is vinyl with finger bumps, and the inside core has foam tubing for gripping a spoon, toothbrush, or other cutlery or utensils.

